

Nest VT™ Vitrification Device

**Secure. Stable. Simple.**



A gloved hand in a dark, textured glove holds a purple pipette. The pipette is angled downwards towards a petri dish containing a white, viscous substance. The background is dark, and the lighting highlights the pipette and the white substance. The overall scene is clinical and precise.

NEST VT VITRIFICATION DEVICE

# A secure solution for embryos and oocytes during vitrification

Securing embryos and oocytes during vitrification begins with you.

The Nest VT Vitrification Device is designed to ease and simplify the process of securing embryos and oocytes during vitrification. With a robust design, including its hallmark “nest” for holding specimens, the Nest VT Vitrification Device helps secure embryos and oocytes for cryo-storage and relocation.

# Nest VT™

## VITRIFICATION DEVICE

---

### Security

- **Recessed “nest” design:** The unique concaved nest for embryos and oocytes helps minimize disturbance to specimens during capping.
- **The confidence of security:** The cap’s friction fit provides tactile feedback so you know specimens are secure inside. The ergonomic, triangular design provides additional stability during the procedure.

### Stability

- **Material integrity:** The device material performs consistently under extreme temperature fluctuations.
- **Purposeful packaging:** The packaging tray keeps the cap and stick separated and protects the device tip from damage during shipping.

### Simplicity

- **Ease of use:** The triangular design provides a stable base and improves device manipulation. A recessed indentation on the cap creates a defined area for forceps, enabling a secure grip during use and capping.
- **Color codes:** Five distinct device colors let you organize and identify specimens the way you prefer.
- **Labeling feature:** The textured and recessed surface area can be used for labeling specimens.






# Ordering information

Order Number	Reference Part Number
G58733	K-NEST-YELLOW
G58734	K-NEST-BLUE
G58735	K-NEST-RED
G58736	K-NEST-GREEN
G58737	K-NEST-PURPLE

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Distribution for details.



# Nest VT™

VITRIFICATION DEVICE

To learn more about the Nest VT Vitrification Device, or to request a demonstration, contact your Cook Medical representative.

### Customer Service

EU Website: [cookmedical.eu](http://cookmedical.eu)  
EDI: [cookmedical.eu/edi](http://cookmedical.eu/edi)  
Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
Belgium: +32 27001702, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
Iceland: +354 800 7615, [is.orders@cookmedical.com](mailto:is.orders@cookmedical.com)  
Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)



USA Website: [cookmedical.com](http://cookmedical.com)  
EDI: [cookmedical.com/edi.do](http://cookmedical.com/edi.do)  
**Americas:**  
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
E-mail: [customersupport@cookmedical.com](mailto:customersupport@cookmedical.com)  
**Australia:**  
Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283  
E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)



AI-ESC-IR-OHNS-PI-RH-SUR-A4