


Procedural accessories for the Zenith Thoracic Endovascular Solution

Needle							
Name	GPN	RPN	Accepts Wire Guide Diameter inch	Gage	Length cm		
One-Part Percutaneous Entry Needle	G04876	SDN-19UT-7.0	.035	19UT	7		
Access Sheaths							
Name	GPN	RPN	Introducer Fr	Accepts Wire Guide Diameter inch	Sheath Length cm	Dilator Length cm	
Performer™ Introducers	G07605	RCF-5.0-35-J	5	.035	13	20	
	G07376	RCF-6.0-35-J	6	.035	13	20	
	G07377	RCF-7.0-35-J	7	.035	13	20	
Introducer Sheaths							
Name	GPN	RPN	Introducer Fr	Accepts Wire Guide Diameter inch	Sheath Length cm	Dilator Length cm	
Large Check-Flo Performer® Introducer with Radiopaque Band	G09691	RCFW-16.0P-38-30-RB	16	.038	30	42	
	G11672	RCFW-18.0P-38-30-RB	18	.038	30	42	
Wire Guides							
Name	GPN	RPN	Diameter inch	Length cm	Flexible Tip Length cm	Tip Curve Radius mm	Coil Length cm
Safe-T-J® Tip Fixed Core Wire Guide	G00902	TSCF-35-180-3	.035	180	2.5	3	N/A
Lunderquist® Extra-Stiff Wire Guide	G31453	TSMG-35-260-LES	.035	260	4	straight	15
	G45352	TSMG-35-300-LES	.035	300	4	straight	15
	G45407	TSCMG-35-300-7-LES	.035	300	4	7.5	15
	G45408	TSCMG-35-300-LESDC	.035	300	4	75/15	15
	G53971	TSCMG-35-300-E-LESDC	.035	300	4	55/15	15
Selective Catheters							
Name	GPN	RPN	Fr	Accepts Wire Guide Diameter inch	Length cm	Tip Configuration	
Beacon® Tip Torcon NB® Advantage Angiographic Catheter	G11457	HNBR5.0-38-100-P-NS-KMP	5	.038	100	KMP	
	G11205	HNBR5.0-38-100-P-NS-MPA	5	.038	100	MPA	
	G11484	HNBR5.0-38-100-P-NS-VERT	5	.038	100	VERT	
Sizing/Flush Catheter							
Name	GPN	RPN	Fr	Accepts Wire Guide Diameter inch	Length cm	Sideports	Number of Bands/ Placement Interval cm
Aurous™ Centimeter Sizing Catheter	G11916	N5.0-35-100-P-10S-PIG-CSC-20	5	.035	100	10	21/1
Dilators							
Name	GPN	RPN	Fr	Accepts Wire Guide Diameter inch	Length cm	Hydrophilic Coating	
Endovascular Dilator Sets	G52833	JCDS-14-EDS-HC	14	.035	45	Yes	
	G52834	JCDS-1618-EDS-HC	16,18	.035	45	Yes	
	G52835	JCDS-2022-EDS-HC	20,22	.035	45	Yes	
	G52836	JCDS-2426-EDS-HC	24,26	.035	45	Yes	
Molding Balloon							
Name	GPN	RPN	Catheter Diameter/ Length Fr/cm	Recommended Introducer Fr	Accepts Wire Guide Diameter inch	Balloon Maximum Inflated Diameter mm	Balloon Maximum Inflation Volume mL
Coda® Balloon Catheter	G23773	CODA-2-10.0-35-140-46	10/140	14	.035	46	60

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 800 7615, is.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com
 EDI: cookmedical.com/edi.do
Americas:
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com
Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com



AI-ESC-IR-OHNS-PI-RH-SUR-A4